

## Needles

	Size	Min. working channel	Max. needle length	Stylet tip shape	Other carachteristics
Cook Echotip ultra ECHO-25	25 G	2.0 mm	80 mm	Round	
Cook Echotip ultra ECHO-1-22	22 G	2.0 mm	80 mm	Round	
Cook Echotip ultra ECHO-3-22	22 G	2.0 mm	80 mm	Round	
Cook Echotip ultra ECHO-19	19 G	2.0 mm	80 mm	Round	
Cook Echotip procore ECHO-HD-25-C	25 G	2.0 mm	80 mm	Round	Biopsy needle
Cook Echotip procore ECHO-HD-22-C	22 G	2.0 mm	80 mm	Round	Biopsy needle
Cook Echotip procore ECHO-HD-19-C	19 G	2.0 mm	80 mm	Round	Biopsy needle
Cook Echotip QuickCore ECHO-HD-19-C	19 G	2.0 mm	80 mm	Beveled	Spring loaded biopsy needle
Cook Echitip ultra ECHO-HD-19 A	19 G	2.0 mm	80 mm	Beveled	Access needle (with atraumatic tip) for EUS-guided drainage
Cook EchoBrush	1 mm	Can be inserted in a Cook Echotip ultra ECHO-19 needle			
Cook Echotip ultra ECHO-20-CPN	20 G	2.8 mm	80 mm	No stylet	Needle with side holes and cone tip for celiac plexus neurolysis
Cook Echotip ECHO-HD-25-EBUS O	25 G	2.0 mm	50 mm	Beveled	EBUS needle for Olympus probes
Cook Echotip ECHO-HD-25-EBUS P	25 G	2.0 mm	50 mm	Beveled	EBUS needle for Pentax probes
Cook Echotip ECHO-HD-22-EBUS O	22 G	2.0 mm	50 mm	Beveled	EBUS needle for Olympus probes
Cook Echotip ECHO-HD-22-EBUS P	22 G	2.0 mm	50 mm	Beveled	EBUS needle for Pentax probes
Cook Echotip procore ECHO-HD-25-EBUS OC	25 G	2.0 mm	50 mm	Beveled	EBUS biopsy needle for Olympus probes

	<b>Size</b>	<b>Min. working channel</b>	<b>Max. needle length</b>	<b>Stylet tip shape</b>	<b>Other carachteristics</b>
Cook Echotip procore ECHO-HD-25-EBUS PC	25 G	2.0 mm	50 mm	Beveled	EBUS biopsy needle for Pentax probes
Cook Echotip procore ECHO-HD-22-EBUS OC	22 G	2.0 mm	50 mm	Beveled	EBUS biopsy needle for Olympus probes
Cook Echotip procore ECHO-HD-22-EBUS PC	22 G	2.0 mm	50 mm	Beveled	EBUS biopsy needle for Pentax probes
Olympus EZ shot NA-200H-8022	22 G	2.8 mm	80 mm	Round	For Olympus probes only
Olympus EZ shot 2 NA-220H-8025	25 G	2.8 mm	80 mm	Round	
Olympus EZ shot 2 NA-220H-8022	22 G	2.8 mm	80 mm	Round	
Olympus EZ shot 2 NA-230H-8022	22 G	2.8 mm	80 mm	Round	Side opening to increase tissue yield
Olympus EZ shot 2 NA-220H-8019	19 G	2.8 mm	80 mm	Round	
Olympus ViziShot NA-201SX-4022	22 G	2.0 mm	40 mm	Round	EBUS needle for Olympus probes
Olympus ViziShot NA-201SX-4021	21 G	2.0 mm	40 mm	Round	EBUS needle for Olympus probes
Boston Scientific Expect	25 G	2.4 mm	80 mm	Beveled	
Boston Scientific Expect	22 G	2.4 mm	80 mm	Beveled	
Boston Scientific Expect	19 G	2.8 mm	80 mm	Beveled	
Boston Scientific Expect Flex	19 G	2.8 mm	80 mm	Beveled	Nitinol needle
Medi Globe SonoTip II	25 G	2.0 mm	85 mm	Round	
Medi Globe SonoTip II	22 G	2.0 mm	85 mm	Round or Beveled	
Medi Globe SonoTip II	19 G	2.8 mm	85 mm	Round	
Medi Globe SonoTip ProControl	25 G	2.0 mm	85 mm	Round	Ergonomic handle
Medi Globe SonoTip	22 G	2.0 mm	85 mm	Round or	Ergonomic handle

	<b>Size</b>	<b>Min. working channel</b>	<b>Max. needle length</b>	<b>Stylet tip shape</b>	<b>Other carachteristics</b>
ProControl				Beveled	
Medi Globe SonoTip ProControl	19 G	2.8 mm	85 mm	Round	Ergonomic handle
Medi Globe SonoTip EBUS Pro	22 G	2.0 mm	40 mm	Round	EBUS needle for Olympus probes
Medi Globe SonoTip EBUS Pro Flex	22 G	2.0 mm	40 mm	Round	EBUS nitinol needle for Olympus probes
Medi Globe SonoTip EBUS Pro Flex	22 G	2.0 mm	40 mm	Round	EBUS nitinol needle for Pentax probes